**NATURE &  
DECOUVERTES**

## **MICROSCOPE ELECTRIQUE**

ELECTRIC MICROSCOPE

MICROSCOPIO ELÉCTRICO

ELEKTRISCHE MICROSCOOP

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**Réf. 42004940**

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**PLEASE CAREFULLY READ THIS MANUAL AND KEEP IT IN A SAFE PLACE**

First, our manufacturer can promise that it is a kind of high-end product with non-polluting and non-toxicity substances, each of whose materials have passed the SGS standards professional certification to be harmless to the human body. The microscope is perfectly suitable for the experimental studies of 8 to 18 years old schoolchildren for the exploration of the microscopic world. It is popular in thousands of preliminary and middle school in Europe and America and is recommended to be the designated product for the "Little Scientist Training Program".


Secondly, we sincerely hope that all the parents and teachers can accompany your children or students to learn the manual together and answer their questions. The warning words are especially provided on the manual as well for the safety and reliability concern of the experiment and machine. Children under 12 years old shall conduct the experiment accompanied by parents or teachers.

At last, we hope that you can use and enjoy our products.

**WARNING**

1. This multifunctional microscope includes some sharp lab wares and some experimental products, so please keep all apparatus well and avoid contact by children.
2. If external power supply is needed for the microscope, please use the product under the supervision of adults.
3. Please do not discard any experimental apparatus to the fire to avoid violent combustion or explosion.
4. It is prohibited to disassemble the machine. In case of any problems, please contact the dealer who will solve the problem together with the manufacturer.
5. It should be noted specially that dry cloth shall be taken to wipe the lens of the machine without using irritant ablent such as soap. The professional cleaning paper or cloth, or the soft rag or paper shall be used to clean the lens of the eyepiece or objective glass.
6. Please put all the items back to the original places after the experiment.
7. In the unpacking and installation process, please pay attention to handle with care without collisions. Each moving parts shall be adjusted property with the entire machine and optical components wiped clean.
8. It should be used rationally and be cleaned after completion of the experiment with the covering of dust cover.
9. The instrument should be placed in dry and clean environment to be away from pollution and corrosive environment.
10. The objective and eyepiece shall be dissembled and placed in dry tanks for long-term non-use. The entire machine shall be covered with the dust cover to avoid getting moldy or contaminated.

## PRODUCT DESCRIPTION

Please look at the picture

- 1 WF10X eyepiece
- 2 Eyepiece sleeve
- 3 Nosepiece
- 4 Objective (4X 10X 40XS)
- 5 Object stage
- 6 Lower light source
- 7 Base
- 8 Curved arm
- 9 Upper light source
- 10 Hand wheel
- 11 Baseplate

## SPECIFICATION

- **Eyepiece:** WF10X
- **Achromatic objective:** 4X, 10X, 40XS
- **Total magnification:**

objective	4X	10X	40X
magnification			
WF10X	40X	100X	400X
- **Observation lens cone:** Monocular, 45° tilt
- **Converter:** three-hole converter

- **Object stage:** 90X90mm metal platform with slice preforming
- **Focusing mechanism:** hand wheel moving
- **Condenser:** NA0.65 Condenser
- **6-hole drive plate light bar:** red, orange, dark green, light green, blue filter.
- **Light source:** LED upper and lower light sources to be adjustable.
- **Power:** 3 AA battery pack (batteries not included)

## OPERATING INSTRUCTION

1. Select the target: the part requiring further observation shall be adjusted to the center with low power lens 4X objective along with the adjustment of the image to the clearest degree for observation at high power lens.
2. The converter shall be turned to replace with the high-powered lens which shall be rotated slowly to observe from the side (to prevent collision with the object slide by high-powered lens), otherwise, it means that the focal length of the low power lens has not been adjusted well. So re-operation shall be conducted.

3. Focal length adjustment: after the converting with high power lens, it should be observed on the eyepiece with the left eye, to see an unclear image, then turn the focusing hand wheel slowly and carefully to achieve a clear image. If the vision brightness is inappropriate, the condenser can be adjusted, and if it is needed to replace the specimen on the slide, the focusing hand wheels wrong direction) to lower the lens stage and remove the specimen on the slide.

## APPLICATION METHOD FOR LOW POWER BIOLOGICAL MICROSCOPE:

**Microscope taking and placing:** the microscope should be stored in a cabinet or case usually to be taken out only in use. The right hand should hold its arm tightly and the left hand should support its base to put the microscope on the experimental stage in front of your left shoulder. It is appropriate that the back-end of the base shall be 1 to 2 cm away from the table edge for the convenience of sitting operation.

**Microscope focusing:** the rotator shall be moved with the thumb and middle (it should not be moved by holding the objective), so that the low power lens can be aligned with the light hole of the low power lens (when the knocking sound is heard in rotation, it means that the objective optical axis has been aligned to the center of the lens cone). Turn on the light source, observe with the left eye in the eyepiece (with the right eye open), while adjusting the dimming knob until the light is even and bright within the view.

**The specimen of the slide shall be placed on the microscope:** a slice specimen shall be put on the stage with the part having cover glass on the upside. It is prohibited to do it in adverse and the part for observation shall be regulated to the right middle side of the light hole.

**Focal length adjustment:** the focusing hand wheel shall be rotated by the left hand in counter-clockwise direction, to lift the platform slowly to about 5mm away from the specimen. It should be noted that observation on the eyepiece shall not be conducted in the platform lifting process. The platform shall be lifted and watched from the right side to avoid over lifting for the damage of the lens or the specimen. Then, both eyes shall be open at the same time with the left eye observing on the eyepiece. The focusing hand wheel shall be rotated slowly by the left-hand clockwise to slow down the platform until the image is clear in vision. In case the image is not in the center of the vision, it can be adjusted (it should be noted that the direction of moving slide it adverse to the moving direction of the image in vision). If the brightness in the vision is inappropriate, the condenser can be lifted and the light rotary knob can be adjusted. If the platform declination has been more than the working distance (>5.40mm) in adjusting the focal length without seeing the image, it means that the operation fails, which shall be conducted for a second time. Be sure not to lift up the microscope stage blindly in a hurry.

## REQUIREMENTS FOR BRIGHTNESS BY THE SPECIMEN SLICE

- When using specimen with good translucency, only the light sources shall be turned on such as (blood slice, plant epidermis slice, etc.).
- When using specimen with poor translucency, both the upper and lower light source shall be turned on such as (root, stem, leaf of plants, chloroplast, fungi, etc.).

## HOW TO MAKE A SPECIMEN

- First the object we prepared shall be very thin so that the light can pass through the specimen and the specimen can be observed by the microscope.
- Sliced onion, leaf epidermis, or saline soaked skin are all specimens easy to be made and are very suitable for beginners.

1. The object slide shall be cleaned.
2. A very thin slice specimen shall be prepared and the ring cutter shall be needed. The preparation process shall be very careful, and you can ask for the help from your parents or teachers.
3. The slice sample shall be clamped by the forceps and put in the center of the slide with a drop of water added by the sucker.
4. The temporary cover glass shall be taken out to be covered to the micro slide gently. Please be careful not to mix any air bubble.
5. The excess water shall be sucked out with an absorbent paper, or it is better to use the paper tower will strong absorbent property.
6. In this way, your temporary sample slice will be completed.

**HOW TO MAKE A PERMANENT SAMPLE SLICE**

1. According to the step 1.2. and 3 for temporary specimen preparation.
2. A few drops of glue or coloring agent shall be added before covering the cover glass.
3. The cover glass shall be covered to the slide gently without air bubble leaving inside.
4. The prepared sample slice shall be placed in safe place after completion for gradually drying.

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EN/WARNING! Not suitable for children under 36 months: risk of swallowing small parts.

WARNING! Discard used batteries immediately. Keep new and used batteries out of the reach of children. If batteries have been swallowed or placed inside any part of the body, seek specialist medical advice immediately. The batteries can only be changed by an adult. Batteries must be disposed of at a dedicated collection point for recycling. The batteries must not be recharged. Do not combine different types of batteries or new and used batteries. Pay attention to the polarity (+ and -) when installing the batteries. Used batteries must be removed. The power terminals (the battery) must not be short-circuited.

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This product must be collected by a specialist service and must not be disposed of with household waste.

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